REG. EU 2017/745 (05/04/2017)  
Norma tecnica EN 14139:2010

MD

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## READING GLASSES

### USE INSTRUCTIONS AND WARNINGS

In order to correctly evaluate the extent and nature of a sight defect, a preliminary sight check is required.

In any case, the evaluation of the sight defect cannot be based of self-diagnosis or on the selection of spectacles considered suitable by following approximations.

The pre-assembled spectacles for presbyopia are constructed industrially with lenses of the same dioptic power, are EXCLUSIVELY recommended for the correction of the "simple" presbyopia.

The pre-assembled spectacles for presbyopia are recommended only for near vision. They must neither be used for driving, nor for far vision or eye protection.

The improper use, i.e. the use of a wrong spectacle gradation, can cause various sight disorders with possible head-ache. In this case, an ophthalmologic visit is required.

The eye specialist only is able to value the state of health of the eyes and the seriousness of visual troubles.

The pre-assembled spectacles for presbyopia must not be used in case of "complex" presbyopia.

The pre-assembled spectacles for presbyopia do not improve sight in case of eye pathologies such as cataract or macular degenerations.

Presbyopia is featured by a physiological and progressive sight weakening that usually reveals after 45 years of age with a series of troubles such as difficulty to read the newspaper or to clearly identify small-size objects from close up.

Presbyopia is considered "simple" when the correction of the sight defect requires the same number of corrective diopters (the diopters are the unit of measure of visual capacity).

Presbyopia is considered "complex" in presence of astigmatism or anisometropia or when the number of diopters of each eye is different. Astigmatism is a sight defect due to an irregular curve of cornea determining an out-of-focus and deformed vision.

Anisometropia is featured by a different refraction defect of the two eyes.

### WARNINGS

- For near vision and reading use only.
- Not for regular use without the approval of an eyecare professional.
- Not for driving or vehicle operation.
- Not for distance vision.
- Not for use as eye protection.
- A regular visit of the eyes from a qualified professional, can determine only your visual necessity and the state of health of your eyes.

Not suitable for driving  
or for use on the road.

The distance between the pupils is of 62 mm.

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